

# ORAL FLUID SWAB

## On-site Testing For Drug Impairment



**Detecting drugs of abuse in oral fluid informs you that the donor is impaired or has recently abused drugs**

A rapid screening test to detect various drugs of abuse in oral fluid.

The Swab range features many innovations to deliver optimal performance and ease-of-use, based on our many years of experience:

- Larger sponge for greater sample volume
- Strainer to extract sample and isolate collector sponge
- Locking lid
- Concave base to deliver sample to test strips
- Multi-line technology to make more efficient use of sample volume
- Hexagonal design makes photocopying or photography easy

### Why Oral Fluid?

When you are workplace or incident screening, you need to know quickly whether someone is under the influence of drugs, or if they have been taken recently.

Oral Fluid is perfect for these instances, because it relates directly to any drugs that are 'live' in blood.

In contrast, someone who has just taken a drug could still record a negative urine test, because the remains of drug take hours to show in a urine sample.

With two components, the Swab is extremely easy and straight forward to use.

Simply place the swab in the mouth, within the teeth to collect the oral fluid, and screw the collector into the tube to activate.

Results are read at 10 minutes, or as soon as a negative result shows.

Two lines indicate a negative result, one line indicates a positive result

### Test Parameters

Drug Type	Retention Time
Amphetamine	72 hrs
Barbiturates	12-50 hrs
Benzodiazepines	24 hrs
Buprenorphine	48-72 hrs
Cannabis	14 hrs
Cocaine	24 hrs
Methadone	48 hrs
Methamphetamines	72 hrs
Opiates	48-72 hrs
Oxycodone	48-72 hrs
Phencyclidine	24 hrs