

INTEGRATED CUP

An all-in-one drug testing solution



Speed and ease of use are of paramount importance when drug screening and the Integrated Cup delivers the most elegant and efficient method of screening.

- One container for collection and testing
- Multipanel Drug Test incorporated
- Easy to read, compact Multiline format
- Tight seal to prevent leaks
- Innovative plunger design secures urine in the cup
- Temperature strip guards against sample adulteration
- Available in many drug combinations

Test up to 12 drug types simultaneously. Available with the following drug groups, and also with Adulteration test parameters:

Drug Type	Code
Amphetamines	AMP
Barbiturates	BAR
Benzodiazepines	BZO
Buprenorphine	BUP
Cannabis	THC
Cocaine	COC
EDDP	EDDP
Fentanyl	FEN
MDMA	MDMA
Methadone	MTD
Methamphetamines	MAMP
Morphine/Opiates	MOP/OPI
Oxycodone	OXY
Phencyclidine	PCP
Propoxyphene	PPX
TCA	TCA
Tramadol	TRA
Oxidants	
Specific Gravity	
pH	
Nitrites	
Glutaraldehyde	
Creatinine	



The hinged lid snaps shut and stores the activation key until you are ready to run the test. The internal rubber gasket gives a tight seal against the lid.



The plunger activates the test and ensures that the correct quantity of sample reaches the test whilst isolating the remainder of the sample. Familiar single drug or multi-line format give an easy to read yes/no result.